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## S.O.P for Reading and Grading Tube Agglutination

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تفسیر واکنش و درجه بندی آگلوتیناسیون Interpretation of Agglutination Reactions	
Macroscopically Observed Findings	Designation
One solid agglutinate	4+
Several large agglutinates	3+
Medium-size agglutinates, clear back ground	2+
Small agglutinates, turbid background	1+
Very small agglutinates, turbid background	1 <sup>w</sup>
Barely visible agglutination, turbid back ground	W <sup>+</sup>
No agglutination	0
Mixtures of agglutinated and unagglutinated red cells (mixed field)	mf
Complete hemolysis	H
Partial hemolysis, some red cells remain	PH

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## References:

1. AABB Technical Manual 16<sup>th</sup> ed. Bethesda, MD 2008
2. Issitt PD, Anstee DJ, Applied blood group serology 4<sup>th</sup> edition, Durham, NC: Montgomery Scientific Press, 1998
3. Judd WJ, Methods in immunohematology 3<sup>rd</sup> ed. Durham, NC: Montgomery Scientific Publication, 2008
4. Denise MH, Modern Blood Banking & Transfusion Practices 5<sup>th</sup> ed. Philadelphia, PA: F.A. Davis company 2005
5. ARC Immunohematology Methods and procedures 1<sup>st</sup> ed. Rockville, MD 1993
6. Guidelines for the Blood Transfusion Services in the UK, 7<sup>th</sup> ed 2005
7. WHO; Safe Blood and Blood products; WHO Distance learning material, WHO PRESS 2004
8. WHO; The Clinical Use of Blood in Medicine, Obstetrics, Pediatrics, Surgery & Anesthesia, Traum & Burns, WHO Press 2002